



Clearly different.
Clearly better.



**NEW! CLIA WAIVED
for Whole Blood!**

Women's Health | HIV / STDs | Respiratory | Cardiovascular | Enterics & GI | General Health

Clearview® COMPLETE HIV 1/2

Clearview® COMPLETE HIV 1/2

A single-use, self-contained, closed system for the collection, processing, and analysis of a whole blood, serum or plasma sample for the detection of HIV 1 and HIV 2 antibodies.



Clearly different

Provides highly sensitive and specific tests for the detection of human immunodeficiency virus (HIV) in whole blood.

- **Rapid Reporting:** Results in 15 minutes.
- **Convenience:** Self-contained system, pouch includes lancet and bandage for fingerstick collection.
- **Reliability:** Built-in control confirms sample addition and proper test performance.
- **Flexibility:** Use with multiple sample types – whole blood (fingerstick or venipuncture), serum or plasma.

Clearly better

- **Reliable results** – 99.7% sensitivity, 99.9% specificity
- **Detects HIV 1 and HIV 2** for a COMPLETE picture in one test.
- **15 minute results** mean STAT decision making, quickly and confidently while patient is present.
- **Minimal sample size of 2.5 µL** reduces sample handling and exposure risk, less invasive for patient.
- **Multiple specimen types** to fit your lab and practice patterns.
- **Long 24-month shelf life** allows your offices to purchase in volume.



Clearview®
COMPLETE HIV 1/2
 25 Test Kit
 Order No. 92111
 CPT Code 86703


Whole Blood – CLIA Waived
 Serum, Plasma –
 Non-waived

Clearview HIV 1/2
Control Kit
 Order No. 92112

CDC Guidelines – Then & Now¹


Then	Now
Targeted testing based on risk assessment	Routine HIV screening recommended for ALL persons ages 13-64 in all healthcare settings Persons at risk should be screened for HIV at least annually
Specific signed consent for HIV testing is required	Opt-out required for screening
HIV screening should be included in the routine prenatal screening tests for all pregnant women	Routine prenatal HIV screening for all pregnant women is recommended, plus repeat 3rd trimester screening for women in high prevalence or high risk areas and pre-labor testing for women with unknown status
Routine counseling (& testing) recommended for high-risk persons and in high prevalence acute-care settings	Prevention counseling encouraged for high risk persons but not required as a part of HIV screening program

Awareness and Serostatus Among People with HIV and Estimates of Transmission²




- People living with HIV/AIDS: 1,039,000 - 1,185,000
- New sexual infections each year: ~32,000

HIV/AIDS Diagnoses among Adults and Adolescents, by Transmission Category 33 States, 2001 – 2004³



Clearview COMPLETE HIV 1/2 Procedure



* Men having sex with men ** Intravenous drug user
 1. CDC. Revised Recommendations for HIV Testing of Adults, Adolescents, and Pregnant Women in Health-Care Settings. MMWR 2006; 55(No.RR-14):1-17
 2. Marks, et al. AIDS 2006; 20:1447-50
 3. MMWR, Nov 18, 2005
 Clearview® is a registered trademark of the Inverness Medical family of companies.

visit
www.invernessmedicalpd.com
 for more information

