



RAPID HIV AND HCV INVALID TEST REPORT
NORTH DAKOTA DEPARTMENT OF HEALTH
DIVISION OF DISEASE CONTROL
SFN 60847 (8/2018)

Counseling, Testing and Referral Site

| | | |
|----------------|------------------|----------|
| Facility Name | Report Date | |
| Street Address | City | ZIP Code |
| Contact Person | Telephone Number | |

Invalid Rapid HIV/HCV Test

| | | |
|----------------------|-----------------------------------------------------------------------------------------------------------|----------------------|
| Date of Invalid Test | Type of Test: <input type="checkbox"/> Chembio SURE CHECK® HIV 1/2 <input type="checkbox"/> OraQuick® HCV | |
| Lot Number | Expiration Date of Test | Date of Last Control |

Invalid Rapid HIV/HCV Test Result

| | |
|-----------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <input type="checkbox"/> Test kit was spilled | <input type="checkbox"/> Test or control line outside valid area (too high or too low) |
| <input type="checkbox"/> Test kit was stored out of temp. range | <input type="checkbox"/> Test or control line did not extend across the window |
| <input type="checkbox"/> Result read too early or too late | <input type="checkbox"/> Other reason (explain below) |
| <input type="checkbox"/> Test kit was expired | <input type="checkbox"/> Reason unknown – please describe appearance of result window (e.g., line at T but not at C, pink result window, etc.). Explain below. |
| <input type="checkbox"/> Forgot to insert a sample | |
| <input type="checkbox"/> No control line | |
| <i>Explain Discrepancy of Test</i> | |
| | |
| | |
| | |
| | |

Quality Assurance

| |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <i>Describe quality assurance follow-up procedures that were conducted to resolve this problem and prevent it from reoccurring (i.e., check storage temperature ranges, inspect expiration dates, and/or verify proper testing techniques).</i> |
| |
| |
| |
| |

Fax this form to 701.328.2499 or mail to:
Attn: HIV.STD.Hepatitis Prevention Coordinator
North Dakota Department of Health
Division of Disease Control
2635 East Main Avenue
Bismarck, N.D. 58506